

ANKOM'S RUMEN SAMPLING BAGS

Standard Fabric:

- *White Polyester Monofilament
- *53 Micron (+/- 10) Pore Size (widely accepted standard)
- *Nitrogen Free

Standard Sizes:

- * 5 x 10 cm bag for 1 gram concentrate sample (10 mg/cm²)
(Actual bag length is 12.5 cm for tying purposes.)
- *10 x 20 cm bag for 5 gram forage sample (12.5 mg/cm²)
(Actual bag length is 22.5 cm for tying purposes.)

Bag Construction:

- *Strong, heat sealed seam that maintains fabric integrity.
- *Bag surface area is designed for commonly used sample sizes.

Prices:

BG505 5x5 cm	\$98.00/100 or \$465.00/500 or \$880.00/1000
BG510 5x10 cm	\$98.00/100 or \$465.00/500 or \$880.00/1000
BG1020 10x20 cm	\$117.00/100 or \$555.00/500 or \$1050.00/1000

The use of a standard rumen bag will eliminate one major variable, thus allowing a more accurate comparison with data obtained at other laboratories using this bag.



www.bardiamond.com & www.bardiamond.info

P.O. Box 60
29865 Sharp Lane
Parma, Idaho 83660-0060 U.S.A.

Phone: 208-722-6761
Fax: 208-722-6686
E-Mail: info@bardiamond.com

Order at www.bardiamond.com