

BAR DIAMOND FLEXIBLE LARGE RUMINANT CANNULAE

Note – STOPPER INCLUDED WITH ALL CANNULAE

See details, prices and place order at www.bardiamond.com

- **Beef, Dairy – all with 4” center diameter; 2-3 lbs each Stopper included.**
 - #1C with rolled inner flange; recommended for dairy breeds \$260
 - #2C with inner flange the same as outer (flat) \$260
 - #3C same as #2C but deeper for thicker (fat) cattle \$270
 - #9C almost impossible to lose but somewhat heavier \$260
 - #9CX same as #9C but for very plump cattle \$280
 - Extra stoppers available for all cannulae \$60-75
 - Add a “U” bolt, or just for fun, add two – inside and out \$20
 - #1W Oblong Washer \$85

- **Beef, Dairy – 3” and 5” center diameter Stopper included.**
 - #4C 3” Use after surgery but before 4” cannula \$250
 - Allows fistula to heal without stretching
 - #1C5 same as #1C but with 5” center diameter; 3 lbs \$290
 - Use when the fistula stretches. Big and heavy.

- **Digestion Bags**
 - BG505 5x5 cm \$98/100; \$465/500 or \$880/1000
 - BG510 5x10 cm \$98/100; \$465/500 or \$880/1000
 - BG1020 10x20 cm \$117/100; \$555/500 or \$1050/1000

- **Other good stuff**
 - #3T Sample Tube to hold digestion bags. Use with 2nd U Bolt \$155
 - #3TP Sample Tube with pump to collect rumen fluid \$255
 - #RT Rumen Fluid Collection Tube attaches to syringe \$150
 - #ST Screw Top added to any stopper. Specify stopper \$85 + stopper
 - #11E-C Esophageal Cannula for cattle \$150
 - #CP1 Hand Soap (you’ll want this!) \$15
 - See web site for intestinal cannulae

- **The Fine Print**

All prices are in US dollars and F.O.B. Parma, Idaho, USA. All orders are shipped by Federal Express unless otherwise requested. Shipping and insurance will be added to the invoice. Buyer assumes all responsibility for use of the products. **March 2022. Prices subject to change without notice.**

FREE T-SHIRT with \$500 rumen cannulae order!

SMALL CANNULAE FOR SHEEP, GOATS, CALVES, LLAMAS, DEER, MUSKOX, WHATEVER!

- **Cannulae 1.5" center diameter. Cannot be removed. ¾ lb Stopper included.**
 - If in doubt, order larger size.
 - **#5C** fits tight. Up to 2" wall thickness. \$240
 - **#6C** up to 3" wall thickness. \$240
 - **#7C** up to 3.5" wall thickness. \$240
 - **Extra stoppers** available for all cannulae \$60

- **Cannulae 2-3" center diameter. Stopper included.**
 - **#8C** 3" center diameter. A small hand will fit through 3" \$240
 - **#10C** 3" rubber. Light weight. Goats. \$185
 - **#11C** 2" rubber. Light weight. Goats. \$165
 - **Extra stoppers** available for all cannulae. \$35-55
 - **Washers** \$85-110

- **Other good stuff**
 - **#11E-S** Esophageal cannula for Sheep \$150
 - **#CP1** Hand Soap (you'll want this!) \$15
 - **See web site for intestinal cannulae**

- **A few comments about small cannulae**
 - Small ruminants are harder to fit than large ruminants.
 - The #8C is best for calves but may be heavy for sheep and is very heavy for goats.
 - Although the #10C and the #11C are very soft and light weight rubber, they do not stay in as well as the plastisol cannulae.
 - Washers are heavy.
 - The #5C, #6C and #7C cannot be removed after inserted because a hand will not fit through the center.
 - See www.bardiamond.com for more detailed information and photos.



www.bardiamond.com & www.bardiamond.info

P.O. Box 60
29865 Sharp Lane
Parma, Idaho 83660-0060 U.S.A.

Phone: 208-722-6761

Fax: 208-722-6686

HELP: info@bardiamond.com